

## REAGENTS AND BUFFERS

### Phosphate-Buffered Saline (PBS) (10X)

80g of NaCl

2.0g of KCl

14.4g of Na<sub>2</sub>HPO<sub>4</sub>

2.4g of KH<sub>2</sub>PO<sub>4</sub>

Mix 800mL ultra-pure water and adjust pH to 7.6 with pure HCl. Top up with ultra-pure water to 1L.

### Phosphate-Buffered Saline w/ Tween20 (PBST)

For 1L: 100mL of PBS 10x + 890mL ultra-pure water + 10 mL Tween 20

### Tris-Buffered Saline (TBS) (10X)

24.23g Trizma HCl

80.06g NaCl

Mix in 800mL ultra-pure water and adjust pH to 7.6 with pure HCl. Top up with ultra-pure water to 1L.

### Tris-Buffered Saline w/ Tween20 (TBST)

For 1L: 100mL of TBS 10x + 890mL ultra-pure water + 10 mL Tween 20

## WESTERN BLOTTING BUFFERS

### RIPA Buffer

150mM NaCl

1.0% NP-40 or 0.1% Triton X-100

0.5% sodium deoxycholate

0.1% SDS (sodium dodecyl sulphate)

50mM Tris-HCl pH8.0

Protease Inhibitors

### Tris-HCl buffer

20mM Tris-HCl pH7.5

Protease Inhibitors

### Denaturing lysis buffer

1% SDS

5mM EDTA

Immediately before use add:

10mM dithiothreitol or beta-mercaptoethanol

Protease inhibitors

15U/mL DNase I

Laemmli 2X buffer/loading buffer

4%SDS  
10% 2-mercaptoethanol  
20% glycerol  
0.004% bromophenol blue  
0.125M Tris-HCl  
Check the pH and adjust pH to 6.8.

Running buffer

25mM Tris base  
190mM glycine  
0.1% SDS  
Check the pH, which should be about pH 8.3. Adjust if necessary.

Transfer buffer (wet)

25mM Tris base  
190mM glycine  
0.1% SDS  
The pH should be about pH 8.3. Adjust if necessary.

Transfer buffer (semi-dry)

48mM Tris  
39mM glycine  
20% methanol  
0.04% SDS

Blocking buffer

5% milk or BSA (bovine serum albumin)  
Add to TBST buffer. Mix well and filter.

## IMMUNOHISTOCHEMISTRY/IMMUNOCYTOCHEMISTRY BUFFERS

Formalin solution (10%)

3.7-4% Formaldehyde (37-40%)  
33mM NaH<sub>2</sub>PO<sub>4</sub>  
46mM Na<sub>2</sub>HPO<sub>4</sub>

Paraformaldehyde (4%)

8% paraformaldehyde

0.2M Phosphate Buffer (PB) pH 7.4

53mM NaH<sub>2</sub>PO<sub>4</sub>  
154mM Na<sub>2</sub>HPO<sub>4</sub>

*Heat 8%PFA solution at 60C while stirring. Once the solution reaches 60C and the PFA dissolves, add 500mL of 0.2M phosphate buffer, to bring the solution to 4%PFA in 0.1M phosphate. Carefully add 1N NaOH until the solution is clear. Cool the solution and filter.*

#### Sodium Citrate Buffer pH6.0

10mM sodium citrate

0.05% Tween20

Mix to dissolve sodium citrate and adjust pH to 6.0 with 1N HCl.

Add Tween20 and mix well.

Store at room temperature for 3 months or at 4C for longer storage.

## ELISA BUFFERS

#### Bicarbonate/carbonate coating buffer (100mM) pH9.6

Antigen or antibody should be diluted in coating buffer to immobilize them to the wells:

29mM Na<sub>2</sub>CO<sub>3</sub>

71mM NaHCO<sub>3</sub>

#### Blocking solution

Commonly used blocking agents are 1% BSA, serum, or non-fat dry milk in PBS.

#### Wash solution

Usually PBS or Tris-buffered saline(pH7.4) with detergent such as 0.05% (v/v) Tween 20 (TBST).

## FLOW CYTOMETRY

#### FACS buffer/antibody dilution buffer

10% FCS

1% sodium azide in PBS

#### Permeabilization

0.1-1% Triton X-100/NP-40 in PBS

#### Fixative

0.01-0.1% paraformaldehyde in PBS.